sonicare

the sonic toothbrush

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
Table of contents

General description
Important safeguards
Dangers
Warnings
Medical warnings
Customer Service Information
Technical Specifications
Preparing for use
Charging your Sonicare®
Using your Sonicare®
Features
Cleaning
Storage
Replacement
Battery recycling
Guarantees & service
Troubleshooting Guide

General description

- | | |
|--------------------------|------------------------------------|
| A Travel cap | F Handle with soft grip |
| B Bristles | G Power on/off button |
| C Brush head | H Dual Speed Control button |
| D Color code ring | I LED charging lights |
| E Removable nut | J Charger |

Please refer to Page 5.

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USE

Sonicare® has been carefully engineered for both effectiveness and safety. However, basic safety precautions should always be followed when using any electrical product, especially when children are present.

DANGERS

To reduce the risk of electrocution or fire:

- 1) DO NOT plug the charger into an incorrect voltage source outlet. It has been designed for use with a specific electrical voltage. Refer to the bottom (or back) of the charger for the correct voltage specification. Voltage converters and plug adapters DO NOT guarantee voltage compatibility. Never force the plug into an outlet.
- 2) DO NOT place or store the charger where it can fall or be pulled into a tub, sink or toilet.
- 3) DO NOT place or drop the charger into water or other liquid.
- 4) DO NOT reach for a charger that has fallen into water. Unplug immediately.

WARNINGS

To reduce the risk of burns, electrocution, fire or injury:

- 1) Close supervision is necessary when this product is used by or on children under the age of 10 or the disabled.
- 2) If you are unable to handle this product properly, discontinue use.
- 3) Use this product only for its intended use as described in this booklet. Do not use attachments other than those recommended by the manufacturer. See voltage information above.
- 4) Never operate the charger if it has a damaged cord or plug, if it has been dropped or if the handle has been damaged in any way. This product contains no serviceable parts. If repair is needed please refer to the Customer Service Information section of this manual and contact one of our representatives.
- 5) Keep the cord away from heated surfaces.
- 6) If your favorite toothpaste includes peroxide, baking soda or other bicarbonate (common in whitening toothpastes), it is important that you thoroughly clean your brush head with soap and water after each use.
- 7) Do not use the charger base outdoors.
- 8) Do not use while bathing.
- 9) Unplug the charger base before cleaning and assure that it is dry before plugging it into an electrical outlet.

MEDICAL WARNINGS

- 1) Consult your dental professional before using this product if you have had oral or gum surgery in the previous 2 months.
- 2) Contact your dental professional if excessive bleeding occurs after using this product or bleeding continues to occur after 1 week of use.
- 3) If you have any questions concerning usage with a pacemaker or other implanted device, consult your physician or the implanted device manufacturer prior to use. Sonicare® has been tested and is compliant with safety standards for electromagnetic devices.
- 4) Avoid using a brush head with crushed or bent bristles. Damaged bristles may break off while brushing. In order to prevent this, replace the brush head with a new one every 6 months or sooner if signs of wear appear.
- 5) The Sonicare® toothbrush is a personal care device and is not intended for use on multiple patients in a dental office or institution.
- 6) This product is designed for cleaning teeth, gums and tongue only. Do not use for any other purpose.

SAVE THESE INSTRUCTIONS

Customer Service Information

If you have any questions after reading this manual, feel free to contact our Customer Service Center or visit our website.

Customer Service Center

1-800-682-7664 (North America), outside North America contact your local Philips Customer Care Center.

0120-418-608 (Japanese language only)

Internet Information

www.sonicare.com (North America) www.philips.com (outside North America)
www.sonicare.jp (Japanese language only)




Technical specifications

Product Name:	Sonicare® Elite e9000 Series
Voltage	100V, 50-60Hz
Speed	Approximately 31,000 brush strokes per minute
Charging Time	24 hours (empty to full)
Operation Time:	A fully charged Sonicare® lasts up to 2 weeks of normal brushing
Battery Type	NiCd
Power Consumption	1.5W max


Preparing for use

To attach the Sonicare® brush head:

- 1** Align the front of the brush head with the front of the toothbrush handle. Match the bars on the brush head with the notches inside the top of the handle and slide the brush head onto the handle.
- 2** Screw the brush head nut clockwise until it is firmly tightened. If the brush head rattles during use, tighten the nut further.
- 3** Remove the travel cap from your brush head when you are ready to brush.
- 4** Use the cord wrap in the bottom of the charger to store excess cord.

Charging your Sonicare®

Sonicare® comes with the batteries completely discharged and should be charged for at least 24 hours during the first few days of use.

▶ It is recommended that you keep your Sonicare® in the charger when not in use to maintain a full battery charge. It takes at least 24 hours to fully charge Sonicare®'s battery.

1 Be sure the bottom of the handle sits back securely in the charger, making complete contact.

▶ Handles of other Sonicare® models do not fit in this charger and will therefore not charge effectively. Do not attempt to force a handle into this charger as this will cause damage.


Using your Sonicare®

Useful tips

- ▶ The brush head comes with a convenient travel cap. Remove the travel cap from the brush head before turning on your Sonicare®.
- ▶ Place the brush head in your mouth before turning on your Sonicare®, and be sure it is turned off before removing it from your mouth, otherwise toothpaste may splatter.
- ▶ Rinse the brush head and bristles after each use.
- ▶ Clean the brush head weekly, including under the nut.
- ▶ A fully charged Sonicare® will provide approximately 2 weeks of typical use (2 minutes of brushing twice a day), making it easy to travel without the charger.


Remember : voltage converters and plug adapters DO NOT guarantee voltage compatibility. Always check for voltage compatibility.

Brushing technique


Follow these simple steps to maximize your brushing experience:

- 1** Before turning the Sonicare® on, place the bristles of the brush head on your gumline at a 45-degree angle.
- 2** Hold the handle with a light grip, turn the Sonicare® on, and gently move the bristles in a small circular motion so the longer bristles reach between your teeth. After a few seconds glide the bristles to the next section. Continue this routine throughout your 2-minute brushing cycle.
- 3** To maximize Sonicare®'s effectiveness apply only light pressure and let the Sonicare® toothbrush do the brushing for you. Do not scrub.

- 4** To make sure you brush evenly throughout the mouth, divide your mouth into 4 sections: outside top, inside top, outside bottom and inside bottom teeth. Begin brushing in section 1 (outside top teeth) and brush for 30 seconds before moving to section 2 (inside top teeth). Continue brushing for 30 seconds in each section progressing through all 4 sections until you've brushed for 2 minutes.
 - 5** After having completed the 2-minute brushing cycle, additional time should be spent brushing the chewing surfaces of your teeth.
 - 6** Press the power button to turn the brush on and off as needed.
- ▶ You may also brush your tongue, with the toothbrush turned on or off, as you prefer.

Tips for specific oral care needs

- ▶ **Stain removal**
In areas where staining occurs, spend extra brushing time to assist with stain removal.
- ▶ **Braces**
The brush should be placed at a 45-degree angle between your braces and gumline. Using light pressure, gently move the bristles in a small circular motion so that the long bristles reach between your teeth for 2 seconds. Then gently turn (roll) the brush handle so the bristles clean the braces and the rest of the tooth before you move the bristles to the next tooth. Continue this until you have brushed all of your teeth. Please note that brush heads will wear out sooner when used on braces.
- ▶ **Dental restorations (dental work)**
Use Sonicare® on dental restorations (such as implants, caps, crowns, bridges and fillings) as you would on your normal tooth surfaces. Work the bristles around all areas of the restoration for thorough cleaning.
- ▶ **Periodontal pockets**
Your dental professional may have indicated that you have periodontal pockets, deep pockets where plaque bacteria collect. If this is the case, after brushing for a full 2 minutes, spend extra time gently brushing these areas to improve your gum health.

Features

Easy-start brushing power ramp-up

- ▶ This Sonicare® model comes with the Easy-start feature activated, which helps you get used to brushing with Sonicare.
- ▶ Brushing power will gently increase to full power over the first 14 brushings.
- ▶ The Easy-start feature is active if you hear 2 beeps immediately after turning on your Sonicare®. Once the appliance has reached optimal brushing power (after the 14th brushing), you will no longer hear the double beep at the beginning of your brushing routine.
- ▶ *Note: Each of the first 14 brushings must be at least 1 minute in length to properly advance through the Easy-start ramp-up cycle.*

To deactivate or activate the Easy-start feature:

- 1** Attach your brush head to the handle.
- 2** Place the handle in the plugged-in charger.
- 3** To deactivate the Easy-start feature: while the handle is still in the charger, press and hold the power button for 5 seconds. You will hear 1 beep to indicate that the Easy-start feature has been deactivated.
- 4** To reactivate the Easy-start feature: while the handle is still in the charger, press and hold the power button for 5 seconds. You will hear 2 beeps to indicate that the Easy-start feature has been activated.


Using the Easy-start feature beyond the initial ramp-up period is not recommended and reduces Sonicare®'s effectiveness in removing plaque.

The Easy-start feature should be deactivated when the Sonicare® is used in clinical trials.

Smartimer 2-minute timer

- ▶ All Sonicare® models come with the Smartimer feature, which automatically sets the brushing time to the 2 minutes recommended by dental professionals.
- ▶ After 2 minutes, your Sonicare® will automatically turn off the toothbrush.
- ▶ If you wish to pause or stop brushing during your 2-minute cycle, press the power button. Press the power button again and the Smartimer will pick up where you left off.
- ▶ *The Smartimer will automatically reset to the beginning of a 2-minute cycle if:*
 - brushing is interrupted for 30 seconds or longer;
 - the handle is placed in the charger.

Quadpacer® 30-second interval timer

- ▶ Quadpacer® makes sure that you brush all sections of your mouth evenly and thoroughly. When using a two-minute brushing cycle, you will hear a short beep and a pause in brushing action at 30, 60 and 90 second intervals.

Programmable Quadpacer® timer


The Programmable Quadpacer® timer offers you the option of customizing the Quadpacer® performance in 1 of 3 ways:

- ▶ You may opt to keep the factory default setting of a 2-minute brushing time with the Quadpacer® feature activated.
- ▶ You may choose to extend your brushing time to 2 1/2 minutes by adding an additional timing interval of 30 seconds with the Quadpacer® feature active.
- ▶ You may choose a 2-minute brushing time, but with the Quadpacer® feature turned off.

To change the Quadpacer® settings:

- 1** Place the handle in the charger while it is plugged into a live electrical outlet.
 - 2** Press the Dual Speed Control button briefly and note which LED lights up.
 - 3** Press and hold the Dual Speed Control button. Continue to hold the button while the LEDs slowly cycle through the timing interval setting options:
- ▶ Green LED in position 4: Quadpacer® active (factory default setting).
 - ▶ Green LED in top (5th) position: 2 1/2 minute setting - Quadpacer® active with additional timing interval of 30 seconds.
 - ▶ Yellow LED in bottom (1st) position: no timing intervals set - Quadpacer® not active.
- 4** To program the desired setting, release the button when the LED corresponding to your choice lights up.

▶ *Note: Selecting 2 1/2 minutes brushing time reduces the number of brushings available from a fully charged battery.*

Dual Speed Control

The Dual Speed Control feature allows you to choose between 2 speed settings: Normal speed (Max) for optimum brushing performance, or Reduced speed (Gentle) when an even more gentle brushing action is desired.

Each time you turn your Sonicare® on, it will begin brushing at the Normal speed setting. To change the speed setting:


- 1** From Normal to Reduced speed: press the Dual Speed Control button once during brushing.
- 2** From Reduced to Normal speed: press the Dual Speed Control button a second time during brushing.

A distinctive two-tone beep indicates that the speed setting has been changed.

It is recommended that you brush at the Normal speed setting to maintain optimal results.

The Normal setting should be selected when the Sonicare® is used in clinical trials.

Battery Charge Indicator

- ▶ If the battery charge of your Sonicare® is low, you will hear 3 beeps after you have brushed your teeth for 2 minutes and the charge indicator LED will blink yellow for 30 seconds.

A low charge indication signifies that 3 or fewer 2-minute brushings remain. When not in use, keep your Sonicare® in the charger to maintain a full battery charge.

Multiple-LED charge indicator

- ▶ When the handle is placed in the charger, the green charge indicator LEDs will blink in a rising pattern as the battery charges. The number of blinking LEDs indicates the level to which the battery has been charged.

When the battery reaches a full charge, the LEDs stop blinking and are illuminated continuously. About once a minute, the charge indicator LEDs go out momentarily and then come back on. This indicates that the battery is being maintained at full charge. The full charge indication will persist until the handle is removed from the charger. The blinking yellow LED indicates that fewer than 3 brushings remain.

Cleaning

By following these simple steps, you should get years of worry-free brushing:

1 Rinse the brush head and bristles after each use, and air dry only.

2 Clean the brush head at least once a week, including under the nut.

3 To remove the nut:

▮ Press down gently on the white area.

Avoid using sharp objects to push on the seal as this may cause damage

▮ Pull the nut over the top of the brush head.


4 To replace the nut, slip the nut over the brush head shaft and pull until it clicks back on.

5 Clean the handle weekly.

6 Clean the charger periodically.

UNPLUG the charger before cleaning. After cleaning the charger, completely dry the charger before plugging it into an electrical outlet.

Mild soap, water and a moist cloth are recommended for cleaning your Sonicare® handle and charger.

Replacement

Brush head

Replace your brush head every 6 months for optimal results. Use only Sonicare® replacement brush heads for this model.

Battery recycling

HANDLE CONTAINS A NICKEL CADMIUM BATTERY.
BATTERY MUST BE RECYCLED OR DISPOSED OF PROPERLY.

- ▶ The battery inside your Sonicare® cannot be replaced, but is easily removed for recycling. At the end of product life and prior to disposal, please remove the battery from the handle by following the instructions below.
- ▶ Do not throw the appliance away with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to preserve the environment.




Removing the batteries

- ▶ Please note that this process is not reversible.

To remove the batteries yourself, you will need a flat-head (standard) screwdriver. Basic safety precautions should be observed when following the procedure outlined below. Be sure to protect your eyes, hands, fingers and the surface on which you are working.

CAUTION

Caution should be exercised when using any household tools (such as a screwdriver) to remove the battery.

- 1** Deplete the battery of any charge by removing the handle from the charger, turning the Sonicare® on and allowing it to run until it stops. Repeat this step until the product will no longer operate.
- 2** Remove the nut from the brush head and place it onto the handle. Tighten the nut on to the handle (A) until the housing comes loose and then pull the nut forward (B).
The handle internals will come out of the housing along with the nut.
- 3** Insert a screwdriver between the electronic circuit board and the plastic part close to each end of the battery; twist and break the connection between the battery and the electronic circuit board in 2 locations.

- 4** Separate the rechargeable battery from the appliance by means of a screwdriver.

US and Canada

For information on proper battery recycling, please visit the Rechargeable Battery Recycling Corporation website at www.rbrc.org or call (678)419-9990. (Philips Oral Healthcare, Inc. is not affiliated with the RBRC.)

GUARANTEES AND SERVICE

SEE SEPARATE PHILIPS WORLDWIDE GUARANTEE LEAFLET FOR STANDARD GUARANTEE.

GUARANTEE RESTRICTIONS

WHAT IS NOT COVERED UNDER WARRANTY:

- BRUSH HEADS
- DAMAGE CAUSED BY MISUSE, ABUSE, NEGLIGENCE OR ALTERATIONS
- NORMAL WEAR AND TEAR, INCLUDING CHIPS, SCRATCHES, ABRASIONS, DISCOLORATION OR FADING

IMPLIED WARRANTIES

ANY IMPLIED WARRANTIES, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED TO THE DURATION OF THE EXPRESS WARRANTIES SET FORTH IN THIS LIMITED WARRANTY. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

LIMITATION OF REMEDIES

IN NO EVENT SHALL PHILIPS ORAL HEALTHCARE, INC. OR ANY OF ITS AFFILIATED OR SUBSIDIARY COMPANIES BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, TORT, OR ANY OTHER LEGAL THEORY. SUCH DAMAGES INCLUDE, WITHOUT LIMITATION, LOSS OF SAVINGS OR REVENUE; LOSS OF PROFIT; LOSS OF USE; THE CLAIMS OF THIRD PARTIES INCLUDING, WITHOUT LIMITATION, DENTISTS AND DENTAL HYGIENISTS; AND COST OF ANY SUBSTITUTE EQUIPMENT OR SERVICES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

Troubleshooting Guide

If you are unable to solve your problem by means of the troubleshooting guide below, in North America contact 1-800-682-7664 or contact us online at www.sonicare.com. Outside North America, contact the Philips Customer Care Center or contact us online at www.philips.com. In Japan, you can contact the Sonicare® Customer Care Center at 0120-418-608 (Japanese language inquiries only)

Problem	Solution
The Sonicare® toothbrush does not work	Recharge your Sonicare® for at least 24 hours.
	Verify your Sonicare® is plugged into a live electrical outlet and not into an outlet controlled by a light switch.
	You may need to re-attach or replace the brush head.
I feel a tickling sensation when I use Sonicare®.	You may experience a slight tickling/tingling sensation when using Sonicare® for the first time. As you become accustomed to brushing with Sonicare®, this sensation will diminish.
The brush head rattles when I use Sonicare®	Sonicare® should produce a gentle hum during operation. Tightening the brush head nut can eliminate excessive rattling. If rattling does not stop you may need to replace the brush head. Make sure you have removed the travel cap before switching the application on.
The brush head chatters against my teeth.	This can be avoided by keeping the brush head properly oriented, with the brush head placed on the gumline at an angle of 45 degrees.
Brushing seems less powerful.	The handle may need recharging.
	The Easy-start feature may be on and needs to be deactivated (see 'Features').
	You may need to clean your brush head.
	You may need to replace your brush head.
	The Dual speed button may be set to reduced speed.
Unpleasant odor	Clean the handle, brush head and charger thoroughly (see 'Cleaning').
	Replace the brush head with a new one.
The LED blinks yellow. What does this mean?	The blinking yellow LED indicates that fewer than 3 brushings remain. Recharge your Sonicare® for at least 24 hours.
The appliance has an extra 30 seconds of brushing time. What is this for?	This function provides additional time for brushing troublesome areas where plaque collects and gum disease starts. It is recommended you use the extra 30 seconds of brushing time to re-brush troublesome areas, massage your gums, or brush your tongue (see 'Features').

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Snoqualmie WA 98065 U.S.A.